


Big Points



Push clear cap in on either side to release the clear cap. Insert the picture. Then replace the clear cap.

How to set up the Big Points:

- Insert a picture into the big points by pushing the clear cap in on either side. This will release the clear cap. Add the picture e.g. picture of Toast. Then replace the clear cap.
- At the back of the device slide the black button onto "Record"



- Press the Big Point down and record your message e.g. Little Mix
- When finished let go of the "Big Point"
- To play your recording slide the black button over to play

When / How to use this:

- Use this throughout the day e.g. to make food choices/music/tv programmes etc. Offer something you know she likes and something she doesn't like.
- Insert the picture under the clear cap – the caps can be removed by pushing in the clear cap at two points on the side
- Replace the clear cap
- Record your message e.g. Little Mix / Classical music. Use single keywords
- Play and show Abigail the two different options

- The option Abigail chooses make sure to then provide e.g. if she picks Classical music, play a short bit of classical music.
- Try to offer the options on either side e.g. left/right. That way you know if Abigail is making a deliberate choice.

